



SOUND REACH

Ultrasound Surgical Equipment
CSUS 8000 Generator and Shears



Tissue Sense
Technology



5mm Sealing



Compact New
Design



Touch Screen
Interface



CE Certified



The Sound Reach™ ultrasonic platform takes a giant leap forward.

The 8000 generator feature a brand new distinctive design and unique features that transforms the way of surgery.



Tissue Sense Technology

A breakthrough technology that changes the way to perform.

Tissue Sense is the combination of both hardware and software allowing the 8000 Generator simultaneously monitor signals from the shears and adjust energy output according to tissue resistance.

This allows a smooth and consistent cutting experience regardless of tissue thickness and default settings, eliminating the need to constantly adjust output settings, so surgeons can now focus on the surgery.



5mm Vessel Sealing

Safe sealing capability for vessels up to 5mm.



Compact Design

A slimmer look to allow more space in the trolley.



Multi-lingual

The 8000 series platform is pre-loaded with multiple languages to support more local users.



Touch Interface

Large touch screen for intuitive and easy user interaction.

Sound Reach Shears

PD Series Shears



Sound Reach Shears offers variety length of shears with multi-purpose design, making it an ideal choice to deliver consistent performance in surgeries.



Curved jaw

Provides great visual clarity and dissection proficiency



Single-hand use

Makes it easy to rotate and activate



Cross-brand Interoperability

Use on your existing generators

Product Code	Shaft Length	Discipline	QTY/BOX
CH14PD	14 cm	Open, ENT procedures	3
CH23PD	23 cm	General, Pediatric procedures	3
CH36PD	36 cm	General, Bariatric, Gynecology, Urology, Colorectal Endoscopic procedures	3
CH45PD	45 cm	Procedures with Obese Patients	3

CSUS 8000 System Components



Generator
CSUS 8000



Transducer
TRA6



Foot switch
FSW2



Trolley (Optional)
CART 6000