

BV-300FR Image Processing Center (Host) Main Technical Parameters

1	Digital Image Processor	High-definition digital image processing center highlights details
2	Image Output	DVI/HDMI
3	White Balance	Automatic white balance
4	Red Gain	-50 - +50 Adjustable
5	Blue Gain	-50 - +50 Adjustable
6	Green Gain	-50 - +50 Adjustable
7	Image brightness	0 - 255 Adjustable
8	Image Saturation	0 - 255 Adjustable
9	Image Contrast	0 - 255 Adjustable
10	Light source type	LED-80W
11	Color Rendering Index	More than 90%
12	Color Temperature	5600K
13	Luminous flux	More than1000lm
14	Voltage	~110V-240V,50Hz/60Hz
15	Air pump pressure	0.06 Mpa
16	airflow	7.5L/min,Suitable for large, medium and small animals
17	Output resolution	HDMI HD video output, 1080P imaging
18	Light source life	20000H
19	Exposure	0-4
20	Tinge	0-360
21	Take photos and videos	High Definition Imaging
22	Memory	Built-in TF card (64G), can record massive videos and take photos, one-click recording and taking photos.
23	Expanded functionality	Purchase a dedicated camera and matching optical mount for hard endoscopes, and the light guide is compatible with laparoscopes and cystoscopes