




Rabbit/Rodent Tabletop Mouth Gag



Speeds up your procedure and save struggling with your operative positioning during all rabbit or rodent dental procedures. Makes examination and treatment much easier and more convenient.





RDT01 Rodent Dental Table Classic

Code RDT01

Specs L 28.2 × W16 cm,
Height of adjustable pole 30 cm

Color Option Yellow, green



RDT021 Rodent Dental Table Multi-functional


Code RDT02

Specs L 26 × W13.2 cm,
Height of adjustable pole 14/17.5 cm

Color Option Blue



Surgical Grade Teflon and Stainless Steel. Specially designed rabbit and rodent table for hands-free dental work on smaller pocket pets and rabbits.



Adjustable platform tilts various angles.



Adjustable bars with gag holds head steady and mouth open for appropriate position without assistance with. Works even if incisors have been removed.



Tootoo Meditech Co., Ltd.

Address: 3A, Phase II, Smart Home,
Baolong Street, Longgang District,
Shenzhen, Guangdong, China

Tel: +86-755-28306385

Email: info@newtootoo.com

Web: www.newtootoo.com

