





VA-900

Anesthesia Ventilator

- Turbine Drive
- Zero Oxygen Consumption
- Up-graded with ICU graded ventilating modes
- 3ml Min Tidal Volume Monitoring
- Widely used Scenarios


PRODUCT STRUCTURE



Pet hospital





Scientific Research



Animal laboratory



Training Institution



PRODUCT FEATURES

- Turbine drive, and no need oxygen source to drive ventilator
- Multiple breathing modes, and ICU graded ventilating modes are safer and healthier for animals
- 3ml Min Tidal Volume Monitoring, provides precise ventilation, more friendly to small animal
- Parameters automatically calculated by only setting the animal's weight
- 9 inch HD touch screen
- More than 20 alarm events, adjustable alarm volume, quick query of alarm records
- Waveforms and rings can be freely matched to make veterinarian's work easier

Parameters

Drive Mode	Turbine and electronic control
Ventilation Mode	VCV, PCV, CPAP/PSV, PRVC, SIMV-PC, VS, SIMV-PRVC, MANUAL
Bellow	0-300ML; 0-1500ML
Waveform	pressure, flow, volume, end expiratory
Loop	pressure-volume/ pressure-flow/flow-volume
Tidal Volume	machine control: 3-1600ml
Manual	5-1600ml
Breathing Rate	1~100bpm
I:E	9.9:1 to 1:9.9
Inhalation Time	0.1~10.0s
Inspiratory Pause%	0-50%
PEEP	OFF, 3~20 cmH ₂ O
Pressure Support	0~60 cmH ₂ O
Pressure Control	5~60 cmH ₂ O
Flow rate Trigger	0.5~20L/min
Pressure Trigger	-1~-20 cmH ₂ O

TooToo Meditech Co., Ltd.

Address: 3A, Phase II, Smart Home, Baolong Street, Longgang District, Shenzhen, Guangdong, China

Tel: +86-755-28306385

Whatsapp: +86 136 3258 4549

Email: info@newtootoo.com

Web: www.newtootoo.com