



Anesthesia Ventilator

VA-200

For Veterinary Use Only

Equipped with 2 kinds of bellows, 300ml and 1500ml. For animals weighing less than 100KG



Real-time display Breath frequency, I:E rate, Tidal volume, Peak airway pressure


Pressure limited working mode to ensure safety

With internal rechargeable battery power supply which can work for over 4 hours normally in full charge.



With AC power failure alarm and battery low voltage alarm

Electronically controlled Mode and Manual mode



VA-200 animal anesthesia ventilator is suitable for automatic respiratory control, respiratory management and respiratory monitoring of different animals during surgery.

Technical Specifications:

Dimension	253mm×224mm×112mm
Weight	3.22kg
Breath frequency	2~120 bpm
I:E ratio	3:1,2:1,1.5:1,1:1,1:1.5,1:2,1:3
Inspiratory tidal volume	20-100-1500ml
Expiratory tidal volume	20~1500ml
Peak inspiratory pressure (PIP) setting	5~60 cmH ₂ O
Airway pressure monitoring	0~99 cmH ₂ O
Animal Weight range	1.5kg-100kg