

TT-VA50

Vertical Autoclave

Suitable for sterilizing of surgical and dental instruments, glassware and biological dressing, food and goods, etc.

Technical Specifications

Parameter	Value
Model	TT-VA50
Capacity	50 L
Working Pressure	0.23 MPa
Working Temp.	105~136 °C
Temp. Precision	0.1 °C
Chamber Material	S304
Chamber Size	φ385 × 515 mm
Noise	≤ 65 dB
Altitude	≤ 2700 m
Standard Accessory	Basket
Consumption	4 kW
External Size (W × D × H)	645 × 610 × 810 mm
Net Weight	91 kg