



TT-SS45

Pulsating Vacuum Pressure Steam Sterilizer

High-capacity table top steam sterilizer with a 45 L chamber and touch screen control panel.



Front view of the unit with the touch screen control panel and the filling port



Specification

Technical parameters of the TT-SS45 model:

Volume Code	TT-SS45
Model (capacity)	45 L
Internal Dimensions	Ø 325 × 595 mm
Equipment Dimensions	775 × 675 × 535 mm
Rated Voltage	AC 110 V / 220 V, 50 ± 1 Hz
Input Power	3000 VA
Sterilization Temperature	121 °C / 134 °C
Circuit Breaker	C16, 400 V~, 50 Hz
Water Tank Capacity	7.5 L (highest water level) 1.0 L (lowest water level)
Noise Limits	< 70 dB (A weighting)
Usage Environment	Ambient temperature: 5 °C – 40 °C Relative humidity: ≤ 85% Atmospheric pressure: 70 kPa – 106 kPa Power supply: a.c. 110 V / 220 V; 50 Hz ± 1 Hz
Net Weight	98 kg
Gross Weight	118 kg
Packaging Dimensions	860 × 710 × 785 mm

Additional Photos



Rear side of the unit: ventilation grille, power inlet and cable outlet



The unit packed in its wooden shipping crate

Contact

TooToo Meditech Co., Ltd.

Address: 3A, Phase II, Smart Home, Baolong Street, Longgang District, Shenzhen, Guangdong, China

Tel: +86-755-28306385

Email: info@newtootoo.com

Web: www.newtootoo.com