

TT-UC20 Ultrasonic Cleaner

Specification (original English text)

Application

Suitable for laboratories, hospital, dentistry, etc. Multifunctional ultrasonic cleaner with LCD display.

Technical data

Parameter	Value
Ultrasonic frequency	40 KHz
Inner liner material	SS304
Capacity	20 L
Adjustable temperature	Room-60 °C
Adjustable time	LCD timer 1-99 minutes
Display	LCD display
Ultrasonic power	360 W
Heating power	500 W
Inner dimensions	330 × 300 × 200 mm (L×W×H)
Outer dimensions	360 × 393 × 445 mm (L×W×H)
N.W.	11.45 kg
G.W.	13 kg

Figure

[FIGURE: Product photograph of the TT-UC20 ultrasonic cleaner.]

The original document contains a product photograph at this position. The photograph shows the unit itself; it carries no labelled controls, connectors or dimension callouts, so it conveys no textual information beyond the technical data above. The image is reproduced below for identification.



Notes on the source text

- The source document prints the ultrasonic power as “36ow”, an evident typing error for 360 W.
- The source document prints the dimension labels as “iner” and “Duter”, evident typing errors for “inner” and “outer”.