

enbio.



Enbio PRO
**Meet the Enbio PRO autoclave.
The Speedmaster.**

Designed
in Switzerland

Another dimension. High performance.

Enbio sterilizers are up to 6 times smaller, 4 times faster and 3 times lighter than a typical class-B autoclave. But when it comes to work efficiency, they keep up the pace.




6x
smaller

4x
faster

3x
lighter

Class B
TripleVacuum

Time is money.



In a busy clinic environment, with dozens of dental or surgical procedures performed daily, the right timing is essential. Enbio changes your sterilization workflow entirely, with reprocessing of expensive instruments between patients, or emergency sterilization available in just minutes.

Enbio PRO

18 min

Typical autoclave

up to 60 min



Enbio PRO. Small dimensions. Big space.

Enbio PRO is an ideal solution as a basic workhorse in a small practice. Superfast sterilization and room to hold most of the instruments, including popular implant cassettes.



Unwrapped FAST / 134°C	Wrapped Class B / 134°C	Wrapped Class B / 121°C	Wrapped / PRION Class B / 134°C
9 min	18 min	30 min	43 min



Unit weight
20 kg

Capacity
Max. packages amount (circa) **16**
Max. recommended batch weight **0.8 kg**
Process chamber volume **5.3 l**

Pressure & temperature
Working pressure **2.1 Bar**
Max. pressure **2.45 Bar**
Max. process tempreature **138°C**

Dimensions
Width **27 cm**
Height **20 cm**
Depth **57.5 cm**


Electrical
Power supply **220-240V / 50Hz**
Nominal power **3.25 kW**
Max. power consumption **16 A**

Misc
Ingress protection **IP20**
Data archiving **USB stick**
Water conductivity* **< 15 µS/cm**

* - demineralized water only


State-of-the-art technology

Enbio's autoclaves are the result of years of extensive research and are a state-of-the-art piece of modern medical technology.



Supereasy.

Plug in. Choose program. Start working.
That's it. No elaborate training.
Just do the job.



1. Plug in
2. Choose program
3. Start working

Installation and operation of Enbio autoclaves
does not require any special training.



A growing network of more than 200 distributors in almost 50 countries works to bring Enbio revolution to your doors. Scan the QR code and find a distributor nearest to you.

Worldwide presence.



Enbio autoclaves are manufactured in the EU and distributed globally through logistics centers covering Europe, North America and the Middle East.

Since 2017, Enbio community has grown to over 60,000 satisfied users who successfully sterilize their tools in our autoclaves.

